

Safety Data Sheet

according to UK REACH Regulation

(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane

Revision date: 23.05.2022

Product code: 4030

Page 1 of 5

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane

Further trade names

8-Azabicyclo[3.2.1]octane-2-carboxylic acid, 3-(4-iodophenyl)-8-methyl, methyl ester, (1R,2S,3S,5S)-

Synonyms:

beta-CIT; 2-CIT

CAS No:

135416-43-2

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**Reference standard for [¹¹C]beta-CITPrecursor for [^{*}]beta-CIT**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds		
	Biomedizinische Forschungsreagenzien GmbH		
Street:	Heinrich-Gläser-Straße 10-14		
Place:	01454 Radeberg		
Telephone:	+49 3528 4041 60	Telefax:	+49 3528 4041 65
e-mail:	info@abx.de		
Contact person:	Dr. Christoph Meyer	Telephone:	+49 3528 4041 8732
e-mail:	meyer@abx.de		
Internet:	http://www.abx.de		

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****GB CLP Regulation**

Acute Tox. 4; H302

Full text of hazard statements: see SECTION 16.

2.2. Label elements**GB CLP Regulation****Hazard components for labelling**

(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane

Signal word: Warning**Pictograms:****Hazard statements**

H302

Harmful if swallowed.

Precautionary statements

P280

Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

Safety Data Sheet

according to UK REACH Regulation

(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane

Revision date: 23.05.2022

Product code: 4030

Page 2 of 5

2.3. Other hazards

Warning - substance not yet tested completely.

SECTION 3: Composition/information on ingredients**3.1. Substances****Chemical characterization**

8-Azabicyclo[3.2.1]octane-2-carboxylic acid, 3-(4-iodophenyl)-8-methyl, methyl ester, (1R,2S,3S,5S)-

Sum formula: C₁₆H₂₀INO₂

Molecular weight: 385.24

Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (GB CLP Regulation)			
135416-43-2	(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane			95 - < 100 %
	Acute Tox. 4; H302			

Full text of H and EUH statements: see section 16.

Specific Conc. Limits, M-factors and ATE

CAS No	EC No	Chemical name	Quantity
		Specific Conc. Limits, M-factors and ATE	
135416-43-2		(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane	95 - < 100 %
		oral: ATE = 500 mg/kg	

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

If swallowed, immediately drink: Water. Call a physician immediately.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**Co-ordinate fire-fighting measures to the fire surroundings. Suitable extinguishing media: Foam. Extinguishing powder. Carbon dioxide (CO₂). Atomized water.**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

Additional information

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures

Safety Data Sheet

according to UK REACH Regulation

(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane

Revision date: 23.05.2022

Product code: 4030

Page 3 of 5

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

Do not breathe dust. Conditions to avoid: skin contact. Eye contact. Wear personal protection equipment.
When using do not eat, drink, smoke, sniff.

Advice on general occupational hygiene

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed.

Further information on storage conditions

Recommended storage temperature:

of °C: 2

up to °C: 8

Protect against: Light.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

8.2. Exposure controls**Appropriate engineering controls**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state:

solid

Colour:

colourless, light yellow

Changes in the physical state

Melting point/freezing point:

110 - 112 °C

Solubility in other solvents

Chloroform.

SECTION 10: Stability and reactivity**10.4. Conditions to avoid**

Light.

10.6. Hazardous decomposition productsCarbon dioxide. Carbon monoxide. Nitrogen oxides (NO_x).hydrogen iodide (HJ).**SECTION 11: Toxicological information**

Safety Data Sheet

according to UK REACH Regulation

(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane

Revision date: 23.05.2022

Product code: 4030

Page 4 of 5

11.1. Information on hazard classes as defined in GB CLP Regulation**Acute toxicity**

CAS No	Chemical name				
	Exposure route	Dose	Species	Source	Method
135416-43-2	(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane				
	oral	ATE 500 mg/kg			

SECTION 12: Ecological information**12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Further information

Do not allow to enter into surface water or drains.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - contaminated packaging

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

Contaminated packaging

Non-contaminated packages may be recycled.

SECTION 14: Transport information**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information

Safety Data Sheet

according to UK REACH Regulation

(-)-2-beta-Carbomethoxy-3-beta-(4-iodophenyl)tropane

Revision date: 23.05.2022

Product code: 4030

Page 5 of 5

Relevant H and EUH statements (number and full text)

H302

Harmful if swallowed.