

Safety Data Sheet

according to UK REACH Regulation

L-Ascorbic acid

Revision date: 17.08.2022

Product code: DP-102-V8

Page 2 of 3

Suitable extinguishing media

Water. Foam. Extinguishing powder. Dry extinguishing powder. Co-ordinate fire-fighting measures to the fire surroundings.

5.3. Advice for firefighters

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General advice

Wear breathing apparatus if exposed to vapours/dusts/aerosols. Suppress gases/vapours/mists with water spray jet.

6.2. Environmental precautions

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on general occupational hygiene

Do not eat, drink, smoke or sneeze at the workplace.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep the packing dry and well sealed to prevent contamination and absorption of humidity.

Further information on storage conditions

Storage up to °C: +30

Protect against: Light.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Additional advice on limit values

To date, no national critical limit values exist.

8.2. Exposure controls

Individual protection measures, such as personal protective equipment

Eye/face protection

Tightly sealed safety glasses.

Hand protection

Single-use gloves. Wash hands before breaks and after work. Protect skin by using skin protective cream.

Skin protection

Wash contaminated clothing prior to re-use.

Respiratory protection

With correct and proper use, and under normal conditions, breathing protection is not required.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state:	solid
Colour:	white
Odour:	odourless

Changes in the physical state

Water solubility: easily soluble.

Safety Data Sheet

according to UK REACH Regulation

L-Ascorbic acid

Revision date: 17.08.2022

Product code: DP-102-V8

Page 3 of 3

Density (at 20 °C):

1,65 g/cm³

SECTION 10: Stability and reactivity

10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide.

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in GB CLP Regulation

SECTION 12: Ecological information

12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations

Dispose of waste according to applicable legislation. Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information

Other applicable information

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

National regulatory information

Water hazard class (D): 1 - slightly hazardous to water

SECTION 16: Other information