

Safety Data Sheet

according to UK REACH Regulation

Tosyl-Fallypride

Revision date: 13.08.2024

Product code: 1550

Page 1 of 4

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Tosyl-Fallypride

Further trade names

Benzamide, 2,3-dimethoxy-5-[3-[[[4-methylphenyl)sulfonyl]oxy]propyl]-N-[[1-(2-propenyl)-2-pyrrolidiny]methyl]-, (S)-

Synonyms:

(S)-2,3-dimethoxy-5-[3-[[[4-methylphenyl)-sulfonyl]oxy]-propyl]-N-[[1-(2-propenyl)-2-pyrrolidiny]methyl]-benzamide

CAS No: 166173-74-6

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Precursor for [¹⁸F]Fallypride

1.3. Details of the supplier of the safety data sheet

Company name:	ABX advanced biochemical compounds Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
E-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
E-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

GB CLP Regulation

This substance is not classified as hazardous in accordance with GB CLP Regulation.

2.2. Label elements

SECTION 3: Composition/information on ingredients

3.1. Substances

Chemical characterization

Benzamide, 2,3-dimethoxy-5-[3-[[[4-methylphenyl)sulfonyl]oxy]propyl]-N-[[1-(2-propenyl)-2-pyrrolidiny]methyl]-, (S)-

Sum formula: C₂₇H₃₆N₂O₆S

Molecular weight: 516.65 g/mol

Relevant ingredients

none (according to UK REACH Regulation)

SECTION 4: First aid measures

4.1. Description of first aid measures

Safety Data Sheet

according to UK REACH Regulation

Tosyl-Fallypride

Revision date: 13.08.2024

Product code: 1550

Page 2 of 4

After inhalation

Provide fresh air.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

If swallowed, immediately drink: Water.

SECTION 5: Firefighting measures

5.1. Extinguishing media**Suitable extinguishing media**Co-ordinate fire-fighting measures to the fire surroundings. Suitable extinguishing media: Foam. Extinguishing powder. Carbon dioxide (CO₂). Atomized water.**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

Additional information

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage

7.1. Precautions for safe handling**Advice on general occupational hygiene**

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container tightly closed.

Further information on storage conditions

Recommended storage temperature:

of °C: -25

up to °C: -15

Store under (Gas): argon. Nitrogen.

Protect against: Light.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters**Additional advice on limit values**

To date, no national critical limit values exist.

8.2. Exposure controls

Safety Data Sheet

according to UK REACH Regulation

Tosyl-Fallypride

Revision date: 13.08.2024

Product code: 1550

Page 3 of 4

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state: liquid viscous
Colour: light yellow
Water solubility: The study does not need to be conducted because the substance is known to be insoluble in water.

Solubility in other solvents

Chloroform Dichloromethane: Diethyl ether ethyl acetate.

SECTION 10: Stability and reactivity

10.4. Conditions to avoid

Light.

10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. Sulfur oxides. Nitrogen oxides (NO_x).

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in GB CLP Regulation

Additional information on tests

The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

SECTION 12: Ecological information

12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Further information

Do not allow to enter into surface water or drains. The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - contaminated packaging

Safety Data Sheet

according to UK REACH Regulation

Tosyl-Fallypride

Revision date: 13.08.2024

Product code: 1550

Page 4 of 4

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

Contaminated packaging

Non-contaminated packages may be recycled.

SECTION 14: Transport information

Other applicable information

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information