

Safety Data Sheet

according to UK REACH Regulation

Bromocholine bromide

Revision date: 31.05.2022

Product code: 6141

Page 1 of 5

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

Bromocholine bromide

Further trade names

(bromomethyl)(2-hydroxyethyl)dimethylammonium bromide

Synonyms:

N,N-Dimethyl-N-(bromomethyl)ethan-2-ol ammonium bromide; N-(Bromomethyl)-2-hydroxy-N,N-dimethylethanaminium bromide; Bromomethylcholine bromide

CAS No: 28508-20-5

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**Reference standard for byproduct of [¹⁸F]Fluoromethylcholine synthesis**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
e-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
e-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****GB CLP Regulation**

Skin Irrit. 2; H315

Eye Irrit. 2; H319

STOT SE 3; H335

Full text of hazard statements: see SECTION 16.

2.2. Label elements**GB CLP Regulation****Hazard components for labelling**

Bromocholine bromide

Signal word: Warning**Pictograms:****Hazard statements**

H315

Causes skin irritation.

H319

Causes serious eye irritation.

H335

May cause respiratory irritation.

Safety Data Sheet

according to UK REACH Regulation

Bromocholine bromide

Revision date: 31.05.2022

Product code: 6141

Page 2 of 5

Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

2.3. Other hazards

Warning - substance not yet tested completely.

SECTION 3: Composition/information on ingredients**3.1. Substances****Chemical characterization**

(bromomethyl)(2-hydroxyethyl)dimethylammonium bromide

Sum formula: C₅H₁₃Br₂NO

Molecular weight: 262.97

Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (GB CLP Regulation)			
28508-20-5	Bromocholine bromide			95 - < 100 %
	Skin Irrit. 2, Eye Irrit. 2, STOT SE 3; H315 H319 H335			

Full text of H and EUH statements: see section 16.

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

If swallowed, immediately drink: Water. Call a physician immediately.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings. Suitable extinguishing media: Foam. Extinguishing powder. Carbon dioxide (CO₂). Atomized water.

5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

Additional information

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

Safety Data Sheet

according to UK REACH Regulation

Bromocholine bromide

Revision date: 31.05.2022

Product code: 6141

Page 3 of 5

6.3. Methods and material for containment and cleaning up**Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

Do not breathe dust. Conditions to avoid: skin contact. Eye contact. Wear personal protection equipment.
When using do not eat, drink, smoke, sniff.

Advice on general occupational hygiene

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed.

Further information on storage conditions

Recommended storage temperature:

of °C: 2

up to °C: 8

Protect against: Light. moisture.

Store under (Gas): argon. Nitrogen.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

8.2. Exposure controls**Appropriate engineering controls**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state:	solid
Colour:	colourless

Changes in the physical state

Melting point/freezing point:	107 - 109 °C
-------------------------------	--------------

Solubility in other solvents

dimethylsulphoxide (DMSO).

SECTION 10: Stability and reactivity**10.4. Conditions to avoid**

Light.

Air, humid.

10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. Nitrogen oxides (NO_x).bromine hydrogen (HBr).

SECTION 11: Toxicological information

Safety Data Sheet

according to UK REACH Regulation

Bromocholine bromide

Revision date: 31.05.2022

Product code: 6141

Page 4 of 5

11.1. Information on hazard classes as defined in GB CLP Regulation**Toxicokinetics, metabolism and distribution**

Toxicological data are not available.

Acute toxicity

Toxicological data are not available.

Further information

Toxicological data are not available.

SECTION 12: Ecological information**12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Further information

Do not allow to enter into surface water or drains.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - contaminated packaging

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

Contaminated packaging

Non-contaminated packages may be recycled.

SECTION 14: Transport information**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information

Safety Data Sheet

according to UK REACH Regulation

Bromocholine bromide

Revision date: 31.05.2022

Product code: 6141

Page 5 of 5

Relevant H and EUH statements (number and full text)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.