

Safety Data Sheet

according to UK REACH Regulation

Sodium hydroxide

Revision date: 07.11.2022

Product code: PEC11-0022-2-V1

Page 1 of 5

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

Sodium hydroxide

Further trade names

Natriumhydroxid

CAS No: 1310-73-2
 Index No: 011-002-00-6
 EC No: 215-185-5

1.2. Relevant identified uses of the substance or mixture and uses advised against**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
e-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
e-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****GB CLP Regulation**

Skin Corr. 1A; H314

Full text of hazard statements: see SECTION 16.

This preparation does not contain harmful substances and/or substances hazardous to the environment in accordance with the substance directive 67/548/EC or substances to which occupational exposure limit values are allocated.

2.2. Label elements**GB CLP Regulation**

Signal word: Danger

Pictograms:

**Hazard statements**

H314 Causes severe skin burns and eye damage.

Precautionary statements

P260 Do not breathe dust/fume/gas/mist/vapours/spray.
 P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
 P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
 P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
 P363 Wash contaminated clothing before reuse.

Safety Data Sheet

according to UK REACH Regulation

Sodium hydroxide

Revision date: 07.11.2022

Product code: PEC11-0022-2-V1

Page 2 of 5

P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P310	Immediately call a POISON CENTER/doctor.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER/doctor.
P405	Store locked up.

SECTION 3: Composition/information on ingredients**3.1. Substances****Hazardous components**

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (GB CLP Regulation)			
1310-73-2	sodium hydroxide; caustic soda			100 %
	215-185-5	011-002-00-6		
	Skin Corr. 1A; H314			

Full text of H and EUH statements: see section 16.

Specific Conc. Limits, M-factors and ATE

CAS No	EC No	Chemical name	Quantity
	Specific Conc. Limits, M-factors and ATE		
1310-73-2	215-185-5	sodium hydroxide; caustic soda	100 %
	Skin Corr. 1A; H314: >= 5 - 100 Skin Corr. 1B; H314: >= 2 - < 5 Skin Irrit. 2; H315: >= 0,5 - < 2 Eye Irrit. 2; H319: >= 0,5 - < 2		

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

If swallowed, immediately drink: Water.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Provide adequate ventilation.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Safety Data Sheet

according to UK REACH Regulation

Sodium hydroxide

Revision date: 07.11.2022

Product code: PEC11-0022-2-V1

Page 3 of 5

Other information

Take up dust-free and set down dust-free.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on general occupational hygiene**

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container tightly closed.

Further information on storage conditions

Maximum storage temperature: +30 °C

Protect against: Light.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Exposure limits (EH40)**

CAS No	Substance	ppm	mg/m ³	fibres/ml	Category	Origin
1310-73-2	Sodium hydroxide	-	2		STEL (15 min)	WEL

8.2. Exposure controls**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state: solid
 Colour: white
 Odour: odourless

Changes in the physical state

Water solubility: easily soluble.

SECTION 10: Stability and reactivity**10.4. Conditions to avoid**

Light.

SECTION 11: Toxicological information**11.1. Information on hazard classes as defined in GB CLP Regulation****SECTION 12: Ecological information****12.1. Toxicity**

CAS No	Chemical name					
	Aquatic toxicity	Dose	[h] [d]	Species	Source	Method
1310-73-2	sodium hydroxide; caustic soda					
	Acute fish toxicity	LC50 mg/l	45,4	96 h	Onchorhynchus mykiss	

12.5. Results of PBT and vPvB assessment

Safety Data Sheet

according to UK REACH Regulation

Sodium hydroxide

Revision date: 07.11.2022

Product code: PEC11-0022-2-V1

Page 4 of 5

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

Contaminated packaging

Water (with cleaning agent). Completely emptied packages can be recycled.

SECTION 14: Transport information**Land transport (ADR/RID)**

14.1. UN number or ID number:	UN 1823
14.2. UN proper shipping name:	Sodium hydroxide, solid
14.3. Transport hazard class(es):	8
14.4. Packing group:	II
Hazard label:	8

Inland waterways transport (ADN)

14.1. UN number or ID number:	UN 1823
14.2. UN proper shipping name:	Sodium hydroxide, solid
14.3. Transport hazard class(es):	8
14.4. Packing group:	II
Hazard label:	8

Marine transport (IMDG)

14.1. UN number or ID number:	UN 1823
14.2. UN proper shipping name:	Sodium hydroxide, solid
14.3. Transport hazard class(es):	8
14.4. Packing group:	II
Hazard label:	8

Air transport (ICAO-TI/IATA-DGR)

14.1. UN number or ID number:	UN 1823
14.2. UN proper shipping name:	Sodium hydroxide, solid
14.3. Transport hazard class(es):	8
14.4. Packing group:	II
Hazard label:	8

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information**

Restrictions on use (REACH, annex XVII):

Entry 3, Entry 75

National regulatory information

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information**Relevant H and EUH statements (number and full text)**

H314 Causes severe skin burns and eye damage.

Safety Data Sheet

according to UK REACH Regulation

Sodium hydroxide

Revision date: 07.11.2022

Product code: PEC11-0022-2-V1

Page 5 of 5

H318

Causes serious eye damage.