

# Safety Data Sheet

according to UK REACH Regulation

**Norverapamil**

Revision date: 31.05.2022

Product code: 5600

Page 1 of 4

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Norverapamil

**Further trade names**

Benzeneacetonitrile, a-[3-[[2-(3,4-dimethoxyphenyl)-ethyl]amino]propyl]-3,4-dimethoxy-a-(1-methylethyl)-

Synonymes:

5-[(3,4-dimethoxyphenethyl)amino]-2-(3,4-dimethoxyphenyl)-2-isopropylvaleronitrile; 2-(3,4-Dimethoxyphenyl)

-2-isopropyl-5-(3,4-dimethoxyphenethylamino)valeronitrile; a-[3-[[2-(3,4-Dimethoxyphenyl)ethyl]amino]propyl]

-3,4-dimethoxy-a-isopropylbenzeneacetonitrile; (±)-Norverapamil; D591

CAS No: 67018-85-3

**1.2. Relevant identified uses of the substance or mixture and uses advised against****Use of the substance/mixture**

Precursor for [11C]Verapamil

**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
e-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
e-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

**1.4. Emergency telephone number:**

+49 3528 4041 60

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****GB CLP Regulation**

This substance is not classified as hazardous in accordance with GB CLP Regulation.

**2.2. Label elements****2.3. Other hazards**

Warning - substance not yet tested completely.

**SECTION 3: Composition/information on ingredients****3.1. Substances****Chemical characterization**

Benzeneacetonitrile, a-[3-[[2-(3,4-dimethoxyphenyl)-ethyl]amino]propyl]-3,4-dimethoxy-a-(1-methylethyl)-

Sum formula: C<sub>26</sub>H<sub>36</sub>N<sub>2</sub>O<sub>4</sub>

Molecular weight: 440.58

**Hazardous components**

none (according to UK REACH Regulation)

**SECTION 4: First aid measures**

**Safety Data Sheet**

according to UK REACH Regulation

**Norverapamil**

Revision date: 31.05.2022

Product code: 5600

Page 2 of 4

**4.1. Description of first aid measures****After inhalation**

Provide fresh air.

**After contact with skin**

After contact with skin, wash immediately with: Water. Change contaminated clothing.

**After contact with eyes**

Rinse immediately carefully and thoroughly with eye-bath or water. In case of troubles or persistent symptoms, consult an ophthalmologist.

**After ingestion**

Rinse mouth immediately and drink plenty of water. Medical treatment necessary.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**Suitable extinguishing media: Water. alcohol resistant foam. Extinguishing powder. Carbon dioxide (CO<sub>2</sub>).**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

**Additional information**

Suppress gases/vapours/mists with water spray jet. Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Wear personal protection equipment.

**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

**6.3. Methods and material for containment and cleaning up****Other information**

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**Advice on general occupational hygiene**

Remove contaminated, saturated clothing immediately. Protect skin by using skin protective cream. After work, wash hands and face. When using do not eat or drink.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Keep container in a well-ventilated place.

**Further information on storage conditions**Recommended storage temperature:  
of °C: -25

**Safety Data Sheet**

according to UK REACH Regulation

**Norverapamil**

Revision date: 31.05.2022

Product code: 5600

Page 3 of 4

up to °C: -15

Protect against: Light.

Store under (Gas): argon. Nitrogen.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

**8.2. Exposure controls****Appropriate engineering controls**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state:

liquid

Colour:

light yellow, light brown

**Changes in the physical state**

Water solubility:

The study does not need to be conducted because the substance is known to be insoluble in water.

**Solubility in other solvents**

Chloroform, Diethyl ether

**SECTION 10: Stability and reactivity****10.4. Conditions to avoid**

Light.

**10.6. Hazardous decomposition products**

Carbon dioxide. Carbon monoxide. nitrogen oxides (NOx).

**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in GB CLP Regulation****Further information**

Toxicological data are not available.

**SECTION 12: Ecological information****12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

**12.6. Endocrine disrupting properties**

This substance does not have endocrine disrupting properties with respect to non-target organisms.

**Further information**

Do not allow to enter into surface water or drains. The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

**SECTION 13: Disposal considerations****13.1. Waste treatment methods**

**Safety Data Sheet**

according to UK REACH Regulation

**Norverapamil**

Revision date: 31.05.2022

Product code: 5600

Page 4 of 4

**Disposal recommendations**

Dispose of waste according to applicable legislation.

**List of Wastes Code - residues/unused products**

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

**List of Wastes Code - used product**

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

**List of Wastes Code - contaminated packaging**

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

**Contaminated packaging**

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

**SECTION 14: Transport information****Other applicable information**

Not a hazardous material with respect to these transportation regulations.

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

**SECTION 16: Other information**