

Safety Data Sheet

according to UK REACH Regulation

N-Trifluoroacetyl-5-acetoxy-2-trimethylstannyl-phenylalanine ethyl ester

Revision date: 10.05.2022

Product code: 3100

Page 1 of 4

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

N-Trifluoroacetyl-5-acetoxy-2-trimethylstannyl-phenylalanine ethyl ester

Further trade names

L-Phenylalanine, 5-(acetoxy)-N-trifluoroacetyl-2-trimethylstannyl, ethyl ester

Synonyms:

N-trifluoroacetyl-3-acetyl-6-trimethylstannyl-L-m-tyrosine ethyl ester;

N-trifluoroacetyl-3-acetoxy-6-trimethylstannyl-L-phenylalanine ethyl ester

CAS No: 148613-10-9

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**

Precursor for [18F]Fluoro-m-L-tyrosine ([18F]FMT)

1.3. Details of the supplier of the safety data sheet

Company name:	ABX advanced biochemical compounds		
	Biomedizinische Forschungsreagenzien GmbH		
Street:	Heinrich-Gläser-Straße 10-14		
Place:	01454 Radeberg		
Telephone:	+49 3528 4041 60	Telefax:	+49 3528 4041 65
e-mail:	info@abx.de		
Contact person:	Dr. Christoph Meyer	Telephone:	+49 3528 4041 8732
e-mail:	meyer@abx.de		
Internet:	http://www.abx.de		

1.4. Emergency telephone number: +49 3528 4041 60**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****GB CLP Regulation**

This substance is not classified as hazardous in accordance with GB CLP Regulation.

2.2. Label elements**2.3. Other hazards**

Warning - substance not yet tested completely.

SECTION 3: Composition/information on ingredients**3.1. Substances****Chemical characterization**

L-Phenylalanine, 5-(acetoxy)-N-trifluoroacetyl-2-trimethylstannyl, ethyl ester

Sum formula: C₁₈H₂₄F₃NO₅Sn

Molecular weight: 510.09

Hazardous components

none (according to UK REACH Regulation)

SECTION 4: First aid measures**4.1. Description of first aid measures**

Safety Data Sheet

according to UK REACH Regulation

N-Trifluoroacetyl-5-acetoxy-2-trimethylstannyl-phenylalanine ethyl ester

Revision date: 10.05.2022

Product code: 3100

Page 2 of 4

After inhalation

Provide fresh air.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

If swallowed, immediately drink: Water.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**Co-ordinate fire-fighting measures to the fire surroundings. Suitable extinguishing media: Foam. Extinguishing powder. Carbon dioxide (CO₂). Atomized water.**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

Additional information

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on general occupational hygiene**

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container tightly closed.

Further information on storage conditions

Recommended storage temperature:

of °C: -25

up to °C: -15

Store under (Gas): argon. Nitrogen.

Protect against: Light.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

8.2. Exposure controls

Safety Data Sheet

according to UK REACH Regulation

N-Trifluoroacetyl-5-acetoxy-2-trimethylstannyl-phenylalanine ethyl ester

Revision date: 10.05.2022

Product code: 3100

Page 3 of 4

Individual protection measures, such as personal protective equipment**Eye/face protection**

Tightly sealed safety glasses.

Hand protection

Single-use gloves. NBR (Nitrile rubber).

Skin protection

Lab apron. Chemical resistant safety shoes.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state: liquid viscous

Colour: light yellow

Changes in the physical state

Water solubility: The study does not need to be conducted because the substance is known to be insoluble in water.

Solubility in other solvents

miscible with most organic solvents

SECTION 10: Stability and reactivity**10.4. Conditions to avoid**

Light.

10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. Nitrogen oxides (NOx). fluorhydric acid. Toxic metal oxide smoke.

SECTION 11: Toxicological information**11.1. Information on hazard classes as defined in GB CLP Regulation****Additional information on tests**

The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

Further information

Toxicological data are not available.

SECTION 12: Ecological information**12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Further information

Do not allow to enter into surface water or drains. The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains.

Safety Data Sheet

according to UK REACH Regulation

N-Trifluoroacetyl-5-acetoxy-2-trimethylstannyl-phenylalanine ethyl ester

Revision date: 10.05.2022

Product code: 3100

Page 4 of 4

List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - contaminated packaging

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

Contaminated packaging

Non-contaminated packages may be recycled.

SECTION 14: Transport information**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information