

Safety Data Sheet

according to UK REACH Regulation

Buffer solution

Revision date: 06.12.2024

Product code: PEST-0132-R-V6

Page 1 of 3

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Buffer solution

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.3. Details of the supplier of the safety data sheet

| | | |
|-----------------|---|-------------------------------|
| Company name: | ABX advanced biochemical compounds Biomedizinische Forschungsreagenzien GmbH | |
| Street: | Heinrich-Gläser-Straße 10-14 | |
| Place: | 01454 Radeberg | |
| Telephone: | +49 3528 4041 60 | Telefax: +49 3528 4041 65 |
| E-mail: | info@abx.de | |
| Contact person: | Dr. Christoph Meyer | Telephone: +49 3528 4041 8732 |
| E-mail: | meyer@abx.de | |
| Internet: | http://www.abx.de | |

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

GB CLP Regulation

This mixture is not classified as hazardous in accordance with GB CLP Regulation.

2.2. Label elements

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Relevant ingredients

none (according to UK REACH Regulation)

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

No special measures are necessary.

After inhalation

Provide fresh air.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

After contact with eyes

If product gets into the eye, keep eyelid open and rinse immediately with large quantities of water, for at least 5 minutes. Subsequently consult an ophthalmologist.

After ingestion

If swallowed, immediately drink: Water.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Safety Data Sheet

according to UK REACH Regulation

Buffer solution

Revision date: 06.12.2024

Product code: PEST-0132-R-V6

Page 2 of 3

Suitable extinguishing media

Water. Foam. Extinguishing powder. Carbon dioxide (CO₂). Co-ordinate fire-fighting measures to the fire surroundings.

5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

Additional information

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used.

Advice on general occupational hygiene

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container tightly closed. Keep locked up. Keep container tightly closed in a cool, well-ventilated place.

Further information on storage conditions

Maximum storage temperature: + 30 °C

Protect against: Light.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Exposure limits (EH40)**

| CAS No | Substance | ppm | mg/m ³ | fibres/ml | Category | Origin |
|-----------|----------------------|-----|-------------------|-----------|---------------|--------|
| 7664-38-2 | Orthophosphoric acid | - | 1 | | TWA (8 h) | WEL |
| | | - | 2 | | STEL (15 min) | WEL |

8.2. Exposure controls**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state: liquid

Colour: colourless

Water solubility: easily soluble.

SECTION 10: Stability and reactivity

Safety Data Sheet

according to UK REACH Regulation

Buffer solution

Revision date: 06.12.2024

Product code: PEST-0132-R-V6

Page 3 of 3

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in GB CLP Regulation

ATEmix calculated

ATE (oral) > 2000 mg/kg; ATE (dermal) > 2000 mg/kg; ATE (inhalation vapour) > 20 mg/l; ATE (inhalation dust/mist) > 5 mg/l

Additional information on tests

The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

SECTION 12: Ecological information

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to UK REACH.

12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

Further information

Do not allow to enter into surface water or drains. The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information

Other applicable information

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

National regulatory information

Water hazard class (D): 1 - slightly hazardous to water

SECTION 16: Other information

Abbreviations and acronyms

Skin Corr: Skin corrosion

(The data for the relevant ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)