

# Safety Data Sheet

according to UK REACH Regulation

## D,L-Homocysteine

Revision date: 12.04.2022

Product code: 2205

Page 1 of 4

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

D,L-Homocysteine

#### Further trade names

Butyric acid, 2-amino-4-mercapto-, DL-

#### Synonyms:

(RS)-Homocysteine; (RS)-2-Amino-4-mercaptobutanoic acid; Butanoic acid, 2-amino-4-mercapto-, (RS)

CAS No: 454-29-5

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.3. Details of the supplier of the safety data sheet

|                 |   |                               |
|-----------------|---|-------------------------------|
| Company name:   | ABX advanced biochemical compounds<br>Biomedizinische Forschungsreagenzien GmbH |                               |
| Street:         | Heinrich-Gläser-Straße 10-14  |                               |
| Place:          | 01454 Radeberg  |                               |
| Telephone:      | +49 3528 4041 60  | Telefax: +49 3528 4041 65     |
| e-mail:         | info@abx.de   |                               |
| Contact person: | Dr. Christoph Meyer   | Telephone: +49 3528 4041 8732 |
| e-mail:         | meyer@abx.de  |                               |
| Internet:       | http://www.abx.de   |                               |

#### 1.4. Emergency telephone number:

+49 3528 4041 60

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### GB CLP Regulation

This substance is not classified as hazardous in accordance with GB CLP Regulation.

#### 2.2. Label elements

##### GB CLP Regulation

##### Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

### SECTION 3: Composition/information on ingredients

#### 3.1. Substances

##### Chemical characterization

Butyric acid, 2-amino-4-mercapto-, DL-

Sum formula: C<sub>4</sub>H<sub>9</sub>NO<sub>2</sub>S

Molecular weight: 135.18

##### Hazardous components

none (according to UK REACH Regulation)

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

# Safety Data Sheet

according to UK REACH Regulation

## D,L-Homocysteine

Revision date: 12.04.2022

Product code: 2205

Page 2 of 4

**After inhalation**

Provide fresh air.

**After contact with skin**

After contact with skin, wash immediately with: Water. Change contaminated clothing.

**After contact with eyes**

Rinse immediately carefully and thoroughly with eye-bath or water.

**After ingestion**

If swallowed, immediately drink: Water. Call a physician immediately.

### SECTION 5: Firefighting measures

**5.1. Extinguishing media****Suitable extinguishing media**Co-ordinate fire-fighting measures to the fire surroundings. Suitable extinguishing media: Foam. Extinguishing powder. Carbon dioxide (CO<sub>2</sub>). Atomized water.**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

**Additional information**

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

### SECTION 6: Accidental release measures

**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

**6.3. Methods and material for containment and cleaning up****Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

### SECTION 7: Handling and storage

**7.1. Precautions for safe handling****Advice on safe handling**

Do not breathe dust. Conditions to avoid: skin contact. Eye contact. Wear personal protection equipment. When using do not eat, drink, smoke, sniff.

**Advice on general occupational hygiene**

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed.

**Further information on storage conditions**

Recommended storage temperature:

of °C: -25

up to °C: -15

Protect against: Light.

### SECTION 8: Exposure controls/personal protection

**8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

# Safety Data Sheet

according to UK REACH Regulation

## D,L-Homocysteine

Revision date: 12.04.2022

Product code: 2205

Page 3 of 4

### 8.2. Exposure controls

#### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state: solid  
Colour: whitish

#### Changes in the physical state

Water solubility: easily soluble.

## SECTION 10: Stability and reactivity

### 10.4. Conditions to avoid

Light.

### 10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. Nitrogen oxides (NO<sub>x</sub>).Sulfur oxides.

## SECTION 11: Toxicological information

### 11.1. Information on hazard classes as defined in GB CLP Regulation

#### Toxicokinetics, metabolism and distribution

Toxicological data are not available.

#### Acute toxicity

Toxicological data are not available.

#### Further information

Toxicological data are not available.

## SECTION 12: Ecological information

### 12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

### 12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

#### Further information

Do not allow to enter into surface water or drains.

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

#### List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

#### List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

## Safety Data Sheet

according to UK REACH Regulation

### D,L-Homocysteine

Revision date: 12.04.2022

Product code: 2205

Page 4 of 4

**List of Wastes Code - contaminated packaging**

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

**Contaminated packaging**

Non-contaminated packages may be recycled.

**SECTION 14: Transport information****Other applicable information**

Not a hazardous material with respect to these transportation regulations.

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

**SECTION 16: Other information**