

# Safety Data Sheet

according to UK REACH Regulation

## n.c.a. Nucleophilic F-L-DOPA precursor (GMP)

Revision date: 03.02.2025

Product code: 1346

Page 1 of 4

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

n.c.a. Nucleophilic F-L-DOPA precursor (GMP)

#### Further trade names

(S)-3-(5-Formyl-4-methoxymethoxy-2-nitro-phenyl)-2-(trityl-amino)-propionic acid tert-butyl ester

#### Synonyms:

ABX128; n.c.a. Nucleophilic F-L-DOPA-Precursor

CAS No: 1614253-58-5

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### Use of the substance/mixture

Precursor for 6-[18F]Fluoro-L-DOPA

#### 1.3. Details of the supplier of the safety data sheet

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
E-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
E-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

#### 1.4. Emergency telephone number:

+49 3528 4041 60

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### GB CLP Regulation

This substance is not classified as hazardous in accordance with GB CLP Regulation.

#### 2.2. Label elements

##### GB CLP Regulation

##### Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

#### 2.3. Other hazards

Warning - substance not yet tested completely.

### SECTION 3: Composition/information on ingredients

#### 3.1. Substances

##### Chemical characterization

(S)-3-(5-Formyl-4-methoxymethoxy-2-nitro-phenyl)-2-(trityl-amino)-propionic acid tert-butyl ester

Sum formula: C35H36N2O7

Molecular weight: 596.67 g/mol

##### Relevant ingredients

none (according to UK REACH Regulation)

# Safety Data Sheet

according to UK REACH Regulation

## n.c.a. Nucleophilic F-L-DOPA precursor (GMP)

Revision date: 03.02.2025

Product code: 1346

Page 2 of 4

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

**After inhalation**

Provide fresh air.

**After contact with skin**

After contact with skin, wash immediately with: Water. Change contaminated clothing.

**After contact with eyes**

Rinse immediately carefully and thoroughly with eye-bath or water.

**After ingestion**

If swallowed, immediately drink: Water. Call a physician immediately.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

**Suitable extinguishing media**Co-ordinate fire-fighting measures to the fire surroundings. Suitable extinguishing media: Foam. Extinguishing powder. Carbon dioxide (CO<sub>2</sub>). Atomized water.

#### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

**Additional information**

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

### SECTION 6: Accidental release measures

#### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

#### 6.3. Methods and material for containment and cleaning up

**Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

**Advice on safe handling**

Do not breathe dust. Conditions to avoid: skin contact. Eye contact. Wear personal protection equipment. When using do not eat, drink, smoke, sniff.

**Advice on general occupational hygiene**

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

#### 7.2. Conditions for safe storage, including any incompatibilities

**Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed.

**Further information on storage conditions**

Recommended storage temperature:

of °C: -25

up to °C: -15

Protect against: Light.

Store under (Gas): argon. Nitrogen.

### SECTION 8: Exposure controls/personal protection

# Safety Data Sheet

according to UK REACH Regulation

## n.c.a. Nucleophilic F-L-DOPA precursor (GMP)

Revision date: 03.02.2025

Product code: 1346

Page 3 of 4

### 8.1. Control parameters

#### Additional advice on limit values

To date, no national critical limit values exist.

### 8.2. Exposure controls

#### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state:	solid
Colour:	light yellow
Solubility in other solvents	
Chloroform, Dichloromethane ,	

## SECTION 10: Stability and reactivity

### 10.4. Conditions to avoid

Light.

### 10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. Nitrogen oxides (NOx).

## SECTION 11: Toxicological information

### 11.1. Information on hazard classes as defined in GB CLP Regulation

### 11.2. Information on other hazards

#### Further information

Toxicological data are not available.

## SECTION 12: Ecological information

### 12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

### 12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

#### Further information

Do not allow to enter into surface water or drains.

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

#### List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

#### List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

## Safety Data Sheet

according to UK REACH Regulation

### n.c.a. Nucleophilic F-L-DOPA precursor (GMP)

Revision date: 03.02.2025

Product code: 1346

Page 4 of 4

**List of Wastes Code - contaminated packaging**

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

**Contaminated packaging**

Non-contaminated packages may be recycled.

**SECTION 14: Transport information****Other applicable information**

Not a hazardous material with respect to these transportation regulations.

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

**SECTION 16: Other information**