

Safety Data Sheet

according to UK REACH Regulation

3'-Fluoro-thymidine (FLT)

Revision date: 29.07.2024

Product code: 1220

Page 1 of 5

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

3'-Fluoro-thymidine (FLT)

Further trade names

Thymidine, 3'-deoxy-3'-fluoro

Synonyms:

3'-Fluoro-3'-deoxythymidine; FLT

Product group: Zwischenprodukt

CAS No: 25526-93-6

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**Reference standard for [18F]FLT
(3'-Deoxy-3'-[18F]fluorothymidine)**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
E-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
E-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

1.4. Emergency telephone number: +49 3528 4041 60**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****GB CLP Regulation**

Acute Tox. 4; H302

Full text of hazard statements: see SECTION 16.

2.2. Label elements**GB CLP Regulation****Signal word:** Warning**Pictograms:****Hazard statements**

H302 Harmful if swallowed.

Precautionary statements

P270	Do not eat, drink or smoke when using this product.
P301+P312	IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
P330	Rinse mouth.

Safety Data Sheet

according to UK REACH Regulation

3'-Fluoro-thymidine (FLT)

Revision date: 29.07.2024

Product code: 1220

Page 2 of 5

SECTION 3: Composition/information on ingredients**3.1. Substances**

Sum formula: C₁₀H₁₃FN₂O₄
 Molecular weight: 244.22 g/mol g/mol

Relevant ingredients

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (GB CLP Regulation)			
25526-93-6	FLT			95 - < 100 %
	Acute Tox. 4; H302			

Full text of H and EUH statements: see section 16.

Specific Conc. Limits, M-factors and ATE

CAS No	EC No	Chemical name	Quantity
	Specific Conc. Limits, M-factors and ATE		
25526-93-6		FLT	95 - < 100 %
	oral: LD50 = 500 mg/kg		

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

Rinse mouth immediately and drink plenty of water. Medical treatment necessary.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**Water spray. Foam. Extinguishing powder. Carbon dioxide (CO₂).**5.2. Special hazards arising from the substance or mixture**In case of fire may be liberated: Carbon dioxide (CO₂). Carbon monoxide Nitrogen oxides (NO_x). Fluorhydric acid.**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

Additional information

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Wear personal protection equipment. Provide adequate ventilation. Avoid dust formation.

Safety Data Sheet

according to UK REACH Regulation

3'-Fluoro-thymidine (FLT)

Revision date: 29.07.2024

Product code: 1220

Page 3 of 5

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Other information

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Do not breathe dust. If handled uncovered, arrangements with local exhaust ventilation have to be used.

Advice on general occupational hygiene

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

Protect skin by using skin protective cream.

Further information on handling

When using do not eat, drink, smoke, sniff.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container dry. Keep container tightly closed.

Further information on storage conditions

storage temperature:

of °C: 2

up to °C: 8

Protect against: Light.

Store under (Gas): argon. Nitrogen.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Additional advice on limit values

To date, no national critical limit values exist.

8.2. Exposure controls

Individual protection measures, such as personal protective equipment

Eye/face protection

Eye protection: Tightly sealed safety glasses.

Hand protection

Tested protective gloves are to be worn: Single-use gloves. NBR (Nitrile rubber).

Skin protection

Body protection: Lab apron. Chemical resistant safety shoes.

Respiratory protection

With correct and proper use, and under normal conditions, breathing protection is not required.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state:

solid

Colour:

colourless

Water solubility:

easily soluble.

Solubility in other solvents

dimethylsulphoxide (DMSO). Methanol.

Safety Data Sheet

according to UK REACH Regulation

3'-Fluoro-thymidine (FLT)

Revision date: 29.07.2024

Product code: 1220

Page 4 of 5

SECTION 10: Stability and reactivity**10.4. Conditions to avoid**

Light.

10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. Fluorines-carbon-hydrogens. Nitrogen oxides (NOx).

SECTION 11: Toxicological information**11.1. Information on hazard classes as defined in GB CLP Regulation****Acute toxicity**

CAS No	Chemical name				
	Exposure route	Dose	Species	Source	Method
25526-93-6	FLT				
	oral	LD50 mg/kg	500	Mouse.	

Carcinogenic/mutagenic/toxic effects for reproduction

No information available.

Additional information on tests

The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

SECTION 12: Ecological information**12.1. Toxicity**

No information available.

12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Further information

Do not allow to enter into surface water or drains. The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - contaminated packaging

Safety Data Sheet

according to UK REACH Regulation

3'-Fluoro-thymidine (FLT)

Revision date: 29.07.2024

Product code: 1220

Page 5 of 5

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information

Other applicable information

No dangerous good in sense of these transport regulations.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information

Abbreviations and acronyms

Acute Tox: Acute toxicity

Relevant H and EUH statements (number and full text)

H302 Harmful if swallowed.