

**Safety Data Sheet**

according to UK REACH Regulation

**Tetrofosmin sulfosalicylate**

Revision date: 20.06.2022

Product code: 7306

Page 1 of 4

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Tetrofosmin sulfosalicylate

**Further trade names**

6,9-bis-(2-ethoxy ethyl)3,12-dioxa-6,9-diphosphatetradecane sulfosalicylate salt (1:2)

Synonymes:

2-[bis(2-ethoxyethyl)phosphanyl]-ethyl-bis(2-ethoxyethyl)phosphane sulfosalicylate salt; Benzoic acid, 2-hydroxy-5-sulfo-, compd. with 6,9-bis(2-ethoxyethyl)-3,12-dioxa-6,9-diphosphatetradecane (2:1)

**1.2. Relevant identified uses of the substance or mixture and uses advised against****Use of the substance/mixture**

Complexation agent for [99mTc]-Labelling

**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds		
	Biomedizinische Forschungsreagenzien GmbH		
Street:	Heinrich-Gläser-Straße 10-14		
Place:	01454 Radeberg		
Telephone:	+49 3528 4041 60	Telefax:	+49 3528 4041 65
e-mail:	info@abx.de		
Contact person:	Dr. Christoph Meyer	Telephone:	+49 3528 4041 8732
e-mail:	meyer@abx.de		
Internet:	http://www.abx.de		
<b>1.4. Emergency telephone number:</b>	+49 3528 4041 60		

**SECTION 2: Hazards identification****2.2. Label elements****GB CLP Regulation****Signal word:** Warning**Pictograms:****Precautionary statements**

P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P280	Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**SECTION 3: Composition/information on ingredients****3.1. Substances**

Sum formula:	C32H52O16P2S2
Molecular weight:	818.83

**Hazardous components**

none (according to UK REACH Regulation)

**Safety Data Sheet**

according to UK REACH Regulation

**Tetrofosmin sulfosalicylate**

Revision date: 20.06.2022

Product code: 7306

Page 2 of 4

**SECTION 4: First aid measures****4.1. Description of first aid measures****After inhalation**

Provide fresh air.

**After contact with skin**

After contact with skin, wash immediately with: Water. Change contaminated clothing.

**After contact with eyes**

Rinse immediately carefully and thoroughly with eye-bath or water. In case of troubles or persistent symptoms, consult an ophthalmologist.

**After ingestion**

Rinse mouth immediately and drink plenty of water. Medical treatment necessary.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**Suitable extinguishing media: Water. alcohol resistant foam. Extinguishing powder. Carbon dioxide (CO<sub>2</sub>).**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

**Additional information**

Suppress gases/vapours/mists with water spray jet. Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Wear personal protection equipment.

**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

**6.3. Methods and material for containment and cleaning up****Other information**

Take up mechanically, placing in appropriate containers for disposal. Treat the recovered material as prescribed in the section on waste disposal.

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**Advice on general occupational hygiene**

Protect skin by using skin protective cream. After work, wash hands and face. When using do not eat or drink.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Keep container in a well-ventilated place.

**Safety Data Sheet**

according to UK REACH Regulation

**Tetrofosmin sulfosalicylate**

Revision date: 20.06.2022

Product code: 7306

Page 3 of 4

**Further information on storage conditions**

Recommended storage temperature:

of °C: 2

up to °C: 8

Protect against: Light.

Store under (Gas): argon. Nitrogen.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

**8.2. Exposure controls****Appropriate engineering controls**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state:

solid

Colour:

colourless

**Changes in the physical state**

Water solubility:

easily soluble.

**Solubility in other solvents**

Water, Ethanol

**SECTION 10: Stability and reactivity****10.4. Conditions to avoid**

Light.

**10.6. Hazardous decomposition products**Sulfur oxides. Carbon dioxide. Carbon monoxide. Nitrogen oxides (NO<sub>x</sub>).**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in GB CLP Regulation****Toxicokinetics, metabolism and distribution**

No information available.

**Acute toxicity**

Acute toxicity, dermal.

**Additional information on tests**

The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

**SECTION 12: Ecological information****12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

**12.6. Endocrine disrupting properties**

This substance does not have endocrine disrupting properties with respect to non-target organisms.

## Safety Data Sheet

according to UK REACH Regulation

### Tetrofosmin sulfosalicylate

Revision date: 20.06.2022

Product code: 7306

Page 4 of 4

#### Further information

Do not allow to enter into surface water or drains. The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Disposal recommendations

Dispose of waste according to applicable legislation.

##### List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

##### List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

##### List of Wastes Code - contaminated packaging

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

##### Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

### SECTION 14: Transport information

#### Land transport (ADR/RID)

##### Other applicable information (land transport)

Not a hazardous material with respect to these transportation regulations.

#### Inland waterways transport (ADN)

##### Other applicable information (inland waterways transport)

Not a hazardous material with respect to these transportation regulations.

#### Marine transport (IMDG)

##### Other applicable information (marine transport)

Not a hazardous material with respect to these transportation regulations.

#### Air transport (ICAO-TI/IATA-DGR)

##### Other applicable information (air transport)

Not a hazardous material with respect to these transportation regulations.

#### Other applicable information

Not a hazardous material with respect to these transportation regulations.

### SECTION 15: Regulatory information

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

##### National regulatory information

Water hazard class (D): 3 - highly hazardous to water

### SECTION 16: Other information