

Safety Data Sheet

according to UK REACH Regulation

Desmethyl-MPT

Revision date: 10.05.2022

Product code: 3260

Page 1 of 4

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

Desmethyl-MPT

Further trade names

1,2,4-Triazine-3,5(2H,4H)-dione, 2-[4-[4-(7-hydroxy-1-naphthalenyl)-1-piperazinyl]butyl]-4-methyl-

Synonyms:

4-Methyl-2-(4-(4-(7-hydroxynaphthalene-1-yl)piperazinyl)butyl)-3,5-dioxo-(2H,4H)-1,2,4-triazine

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**

Precursor for [11C]MPT

1.3. Details of the supplier of the safety data sheet

Company name:	ABX advanced biochemical compounds		
	Biomedizinische Forschungsreagenzien GmbH		
Street:	Heinrich-Gläser-Straße 10-14		
Place:	01454 Radeberg		
Telephone:	+49 3528 4041 60	Telefax:	+49 3528 4041 65
e-mail:	info@abx.de		
Contact person:	Dr. Christoph Meyer	Telephone:	+49 3528 4041 8732
e-mail:	meyer@abx.de		
Internet:	http://www.abx.de		

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****GB CLP Regulation**

This substance is not classified as hazardous in accordance with GB CLP Regulation.

2.2. Label elements**GB CLP Regulation****Special labelling of certain mixtures**

Caution - substance not yet tested completely.

SECTION 3: Composition/information on ingredients**3.1. Substances****Chemical characterization**

1,2,4-Triazine-3,5(2H,4H)-dione, 2-[4-[4-(7-hydroxy-1-naphthalenyl)-1-piperazinyl]butyl]-4-methyl-

Sum formula: C₂₂H₂₇N₅O₃

Molecular weight: 409.48

Hazardous components

none (according to UK REACH Regulation)

SECTION 4: First aid measures**4.1. Description of first aid measures**

Safety Data Sheet

according to UK REACH Regulation

Desmethyl-MPT

Revision date: 10.05.2022

Product code: 3260

Page 2 of 4

After inhalation

In case of breathing difficulties administer oxygen. Provide fresh air. In all cases of doubt, or when symptoms persist, seek medical advice.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing. Medical treatment necessary.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of troubles or persistent symptoms, consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink plenty of water. Medical treatment necessary.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Suitable extinguishing media: Water. Foam. Extinguishing powder. Carbon dioxide (CO₂).

5.2. Special hazards arising from the substance or mixture

In case of fire may be liberated: Nitrogen oxides (NO_x). Carbon monoxide Carbon dioxide (CO₂).

Additional information

Suppress gases/vapours/mists with water spray jet. Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Provide adequate ventilation. Avoid dust formation. Do not breathe dust. Avoid contact with skin, eyes and clothes. Wear personal protection equipment.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Avoid generation of dust. Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Protect skin by using skin protective cream. After work, wash hands and face. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Keep container in a well-ventilated place.

Further information on storage conditions

Recommended storage temperature: :
of °C: -25 up to °C: -15

Safety Data Sheet

according to UK REACH Regulation

Desmethyl-MPT

Revision date: 10.05.2022

Product code: 3260

Page 3 of 4

Protect against: Light.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

8.2. Exposure controls**Appropriate engineering controls**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

Individual protection measures, such as personal protective equipment**Eye/face protection**

Eye protection: Tightly sealed safety glasses.

Hand protection

Tested protective gloves are to be worn: Single-use gloves. NBR (Nitrile rubber).

Skin protection

Body protection: Lab apron. Chemical resistant safety shoes.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state: solid

Colour: whitish

Changes in the physical state**Solubility in other solvents**

Methanol. dimethylsulphoxide (DMSO).

SECTION 10: Stability and reactivity**10.4. Conditions to avoid**

Light.

SECTION 11: Toxicological information**11.1. Information on hazard classes as defined in GB CLP Regulation****Acute toxicity**

Acute toxicity, oral. Acute toxicity, inhalant. Acute toxicity, dermal.

Additional information on tests

The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

Further information

Toxicological data are not available.

SECTION 12: Ecological information**12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Further information

Do not allow to enter into surface water or drains. The classification was carried out according to the calculation method of the Preparations Directive (1999/45/EC).

Safety Data Sheet

according to UK REACH Regulation

Desmethyl-MPT

Revision date: 10.05.2022

Product code: 3260

Page 4 of 4

SECTION 13: Disposal considerations**13.1. Waste treatment methods****List of Wastes Code - residues/unused products**

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

List of Wastes Code - contaminated packaging

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information