

Safety Data Sheet

according to UK REACH Regulation

TET

Revision date: 09.05.2022

Product code: 3050

Page 1 of 4

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

TET

Further trade names

L-Tyrosine, O-[2-[[[(4-methylphenyl)sulfonyl]oxy]ethyl]-N-(triphenylmethyl)-1,1-dimethylethyl ester

Synonyme: L-Tyrosine, O-(2-tosyloxyethyl)-N-trityl, tert-butyl ester; (2S)-O-(2'-tosyloxyethyl)-N-trityl-tyrosine-tert-butyl ester

CAS No: 478037-15-9

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**Precursor für O-(2-[¹⁸F]Fluoroethyl)-L-tyrosine**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
e-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
e-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

1.4. Emergency telephone number: +49 3528 4041 60**SECTION 2: Hazards identification****2.2. Label elements****GB CLP Regulation****Hazard statements**

H302+H312+H332 Harmful if swallowed, in contact with skin or if inhaled.

Precautionary statementsP280 Wear eye protection/face protection.
P280 Wear protective gloves.**2.3. Other hazards**

Warning - substance not yet tested completely.

SECTION 3: Composition/information on ingredients**3.1. Substances**Sum formula: C₄₁H₄₃NO₆S

Molecular weight: 677.85

Hazardous components

none (according to UK REACH Regulation)

SECTION 4: First aid measures**4.1. Description of first aid measures**

Safety Data Sheet

according to UK REACH Regulation

TET

Revision date: 09.05.2022

Product code: 3050

Page 2 of 4

After inhalation

In case of breathing difficulties administer oxygen. Provide fresh air. In all cases of doubt, or when symptoms persist, seek medical advice.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing. Medical treatment necessary.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of troubles or persistent symptoms, consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink plenty of water. Medical treatment necessary.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Suitable extinguishing media: Water. Foam. Extinguishing powder. Carbon dioxide (CO₂).

5.2. Special hazards arising from the substance or mixture

In case of fire may be liberated: Nitrogen oxides (NO_x). Carbon monoxide Carbon dioxide (CO₂). Sulfur oxides.

Additional information

Suppress gases/vapours/mists with water spray jet. Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Provide adequate ventilation. Avoid dust formation. Do not breathe dust. Avoid contact with skin, eyes and clothes. Wear personal protection equipment.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Avoid generation of dust. Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Protect skin by using skin protective cream. After work, wash hands and face. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Keep container in a well-ventilated place.

Further information on storage conditions

Recommended storage temperature: :

Safety Data Sheet

according to UK REACH Regulation

TET

Revision date: 09.05.2022

Product code: 3050

Page 3 of 4

of °C: -25 up to °C: -15

Protect against: Light.

Store under (Gas): argon. Nitrogen.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Additional advice on limit values

To date, no national critical limit values exist.

8.2. Exposure controls

Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

Individual protection measures, such as personal protective equipment

Eye/face protection

Eye protection: Tightly sealed safety glasses.

Hand protection

Tested protective gloves are to be worn: Single-use gloves. NBR (Nitrile rubber).

Skin protection

Body protection: Lab apron. Chemical resistant safety shoes.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state: solid
Colour: whitish

Changes in the physical state

Water solubility: The study does not need to be conducted because the substance is known to be insoluble in water.

SECTION 10: Stability and reactivity

10.4. Conditions to avoid

Light.

10.5. Incompatible materials

Water. Oxidizing agents, strong.

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in GB CLP Regulation

Acute toxicity

Acute toxicity, oral. Acute toxicity, inhalant. Acute toxicity, dermal.

Further information

Toxicological data are not available.

SECTION 12: Ecological information

12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Safety Data Sheet

according to UK REACH Regulation

TET

Revision date: 09.05.2022

Product code: 3050

Page 4 of 4

Further information

Do not allow to enter into surface water or drains.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Dispose of waste according to applicable legislation.

Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information**Relevant H and EUH statements (number and full text)**

H302+H312+H332 Harmful if swallowed, in contact with skin or if inhaled.