

# Safety Data Sheet

according to UK REACH Regulation

## Dimethyl-8-acetyl-hexahydropyrrolo[2,3]indol-1,2-dicarboxylate

Revision date: 01.06.2022

Product code: 6500

Page 1 of 4

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Dimethyl-8-acetyl-hexahydropyrrolo[2,3]indol-1,2-dicarboxylate

#### Further trade names

Pyrrolo[2,3-b]indole-1,2(2H)-dicarboxylic acid, 8-acetyl-3,3a,8,8a-tetrahydro-, dimethyl ester, (2a,3a beta,8a beta)-

#### Synonyms:

Dimethyl-(2S, 3aR, 8aS)-8-acetyl-1,2,3,3a,8a-hexahydropyrrolo[2,3]indol-1,2-dicarboxylate

CAS No: 79465-83-1

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### Use of the substance/mixture

Precursor for alpha-[11C]Methyl-L-tryptophan

#### 1.3. Details of the supplier of the safety data sheet

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
e-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
e-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

#### 1.4. Emergency telephone number:

+49 3528 4041 60

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### GB CLP Regulation

This substance is not classified as hazardous in accordance with GB CLP Regulation.

#### 2.2. Label elements

##### GB CLP Regulation

##### Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

#### 2.3. Other hazards

Warning - substance not yet tested completely.

### SECTION 3: Composition/information on ingredients

#### 3.1. Substances

##### Chemical characterization

Pyrrolo[2,3-b]indole-1,2(2H)-dicarboxylic acid, 8-acetyl-3,3a,8,8a-tetrahydro-, dimethyl ester, (2a,3a beta,8a beta)-

Sum formula: C<sub>16</sub>H<sub>18</sub>N<sub>2</sub>O<sub>5</sub>

Molecular weight: 318.32

##### Hazardous components

none (according to UK REACH Regulation)

# Safety Data Sheet

according to UK REACH Regulation

## Dimethyl-8-acetyl-hexahydropyrrolo[2,3]indol-1,2-dicarboxylate

Revision date: 01.06.2022

Product code: 6500

Page 2 of 4

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

**After inhalation**

Provide fresh air.

**After contact with skin**

After contact with skin, wash immediately with: Water. Change contaminated clothing.

**After contact with eyes**

Rinse immediately carefully and thoroughly with eye-bath or water.

**After ingestion**

If swallowed, immediately drink: Water. Call a physician immediately.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

**Suitable extinguishing media**Co-ordinate fire-fighting measures to the fire surroundings. Suitable extinguishing media: Foam. Extinguishing powder. Carbon dioxide (CO<sub>2</sub>). Atomized water.

#### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

**Additional information**

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

### SECTION 6: Accidental release measures

#### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

#### 6.3. Methods and material for containment and cleaning up

**Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

**Advice on safe handling**

Do not breathe dust. Conditions to avoid: skin contact. Eye contact. Wear personal protection equipment. When using do not eat, drink, smoke, sniff.

**Advice on general occupational hygiene**

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

#### 7.2. Conditions for safe storage, including any incompatibilities

**Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed.

**Further information on storage conditions**

Recommended storage temperature:

of °C: 2

up to °C: 8

Protect against: Light.

Store under (Gas): argon. Nitrogen.

# Safety Data Sheet

according to UK REACH Regulation

## Dimethyl-8-acetyl-hexahydropyrrolo[2,3]indol-1,2-dicarboxylate

Revision date: 01.06.2022

Product code: 6500

Page 3 of 4

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

##### Additional advice on limit values

To date, no national critical limit values exist.

#### 8.2. Exposure controls

##### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state: solid  
Colour: colourless

##### Changes in the physical state

Melting point/freezing point: 148 - 152 °C

##### Solubility in other solvents

Methanol. Dichloromethane: Acetone

### SECTION 10: Stability and reactivity

#### 10.4. Conditions to avoid

Light.

#### 10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. Nitrogen oxides (NO<sub>x</sub>).

### SECTION 11: Toxicological information

#### 11.1. Information on hazard classes as defined in GB CLP Regulation

##### Information on likely routes of exposure

Toxicological data are not available.

##### Further information

Toxicological data are not available.

### SECTION 12: Ecological information

#### 12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

#### 12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

##### Further information

Do not allow to enter into surface water or drains.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

##### List of Wastes Code - residues/unused products

## Safety Data Sheet

according to UK REACH Regulation

### Dimethyl-8-acetyl-hexahydropyrrolo[2,3]indol-1,2-dicarboxylate

Revision date: 01.06.2022

Product code: 6500

Page 4 of 4

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

**List of Wastes Code - used product**

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

**List of Wastes Code - contaminated packaging**

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

**Contaminated packaging**

Non-contaminated packages may be recycled.

#### SECTION 14: Transport information

**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

#### SECTION 15: Regulatory information

**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

#### SECTION 16: Other information