

Safety Data Sheet

according to UK REACH Regulation

FBnTP Precursor

Revision date: 22.04.2025

Product code: 7058

Page 1 of 6

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

FBnTP Precursor

Further trade names

Triphenyl(4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)benzyl)phosphonium trifluoromethanesulfonate

synonym: Phosphonium, triphenyl[[4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)phenyl]methyl]-, 1,1,1-trifluoromethanesulfonate (1:1) (ACI)

Substance name: FBnTP Precursor

CAS No: 2000202-05-9

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**Precursor for [¹⁸F]Fluorobenzyltriphenylphosphonium cation**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
E-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
E-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****GB CLP Regulation**

Skin Irrit. 2; H315

Eye Irrit. 2; H319

Full text of hazard statements: see SECTION 16.

2.2. Label elements**GB CLP Regulation****Signal word:**

Warning

Pictograms:**Hazard statements**

H315

Causes skin irritation.

H319

Causes serious eye irritation.

Precautionary statements

P264

Wash hands thoroughly after handling.

P280

Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

Safety Data Sheet

according to UK REACH Regulation

FBnTP Precursor

Revision date: 22.04.2025

Product code: 7058

Page 2 of 6

P302+P352	IF ON SKIN: Wash with plenty of water and soap.
P321	Specific treatment (see ... on this label).
P332+P313	If skin irritation occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.

2.3. Other hazards

This substance does not meet the criteria for classification as PBT or vPvB.

SECTION 3: Composition/information on ingredients**3.1. Substances****Chemical characterization**

Triphenyl(4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)benzyl)phosphonium trifluoromethanesulfonate

Sum formula: C₃₂H₃₃BF₃O₅PS

Molecular weight: 628.45 g/mol

Relevant ingredients

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (GB CLP Regulation)			
2000202-05-9	FBnTP Precursor			95 - < 100 %
	Skin Irrit. 2, Eye Irrit. 2; H315 H319			

Full text of H and EUH statements: see section 16.

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air.

If breathing is irregular or stopped, administer artificial respiration.

When in doubt or if symptoms are observed, get medical advice.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

In case of skin irritation, seek medical treatment.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

If eye irritation persists: Get medical advice/attention.

After ingestion

Rinse mouth immediately and drink plenty of water.

IF SWALLOWED: Call a doctor if you feel unwell.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**Water spray., Extinguishing powder., Carbon dioxide (CO₂)., alcohol resistant foam.

Safety Data Sheet

according to UK REACH Regulation

FBnTP Precursor

Revision date: 22.04.2025

Product code: 7058

Page 3 of 6

5.2. Special hazards arising from the substance or mixture

Thermal decomposition can lead to the escape of irritating gases and vapours.

5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Wear personal protection equipment. Provide adequate ventilation. Avoid dust formation.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal. , Keep container tightly closed.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

Do not breathe dust. If handled uncovered, arrangements with local exhaust ventilation have to be used.

Avoid contact with skin, eyes and clothes.

Avoid breathing dust/mist.

Advice on general occupational hygiene

Change contaminated clothing.

Wash hands before breaks and after work.

When using do not eat or drink.

Protect skin by using skin protective cream.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container dry.

Store in a well-ventilated place. Keep container tightly closed.

Further information on storage conditions

storage temperature:

of °C: -25

up to °C: -15

Protect against: Light.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Additional advice on limit values**

This mixture does not contain any substances which either present a health or environmental hazard according to Regulation (EC) No. 1272/2008 or have an occupational exposure limit assigned.

8.2. Exposure controls**Appropriate engineering controls**

Technical ventilation of workplace, If handled uncovered, arrangements with local exhaust ventilation have to be used.

Individual protection measures, such as personal protective equipment

Safety Data Sheet

according to UK REACH Regulation

FBnTP Precursor

Revision date: 22.04.2025

Product code: 7058

Page 4 of 6

Eye/face protection

Eye protection: Tightly sealed safety glasses.

Hand protection

Tested protective gloves are to be worn: Single-use gloves.

Skin protection

Body protection: Lab apron. Chemical resistant safety shoes.

Respiratory protection

Usually no personal respirative protection necessary.

If handled uncovered, arrangements with local exhaust ventilation have to be used.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state: solid
Colour: colourless , beige
Solubility in other solvents
Dichloromethane, dimethylsulphoxide (DMSO). , Dimethylformamide

SECTION 10: Stability and reactivity**10.2. Chemical stability**

The product is chemically stable under recommended conditions of storage, use and temperature.

10.4. Conditions to avoid

Heat, Keep away from sources of heat (e.g. hot surfaces), sparks and open flames. , Protect from sunlight.

10.5. Incompatible materials

Oxidizing agents, strong. , water

10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. Nitrogen oxides (NOx).

SECTION 11: Toxicological information**11.1. Information on hazard classes as defined in GB CLP Regulation****Toxicokinetics, metabolism and distribution**

No information available.

Acute toxicity

Based on available data, the classification criteria are not met.

Irritation and corrosivity

Skin corrosion/irritation: Causes skin irritation.

Serious eye damage/eye irritation: Causes serious eye irritation.

Sensitising effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Germ cell mutagenicity: Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

No information available.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Safety Data Sheet

according to UK REACH Regulation

FBnTP Precursor

Revision date: 22.04.2025

Product code: 7058

Page 5 of 6

Aspiration hazard

Based on available data, the classification criteria are not met.

11.2. Information on other hazards**Further information**

Toxicological data are not available.

SECTION 12: Ecological information**12.1. Toxicity**

Based on available data, the classification criteria are not met.

No information available.

12.2. Persistence and degradability

No information available.

12.3. Bioaccumulative potential

No information available.

12.4. Mobility in soil

No information available.

12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

No information available.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

No information available.

12.7. Other adverse effects

No information available.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

Hazardous waste according to Directive 2008/98/EC (waste framework directive).

The waste code has to be identified in agreement with the disposal company or the competent authority.

List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - contaminated packaging

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information

Safety Data Sheet

according to UK REACH Regulation

FBnTP Precursor

Revision date: 22.04.2025

Product code: 7058

Page 6 of 6

Other applicable information

No dangerous good in sense of these transport regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

15.2. Chemical safety assessment

For this substance a chemical safety assessment has not been carried out.

SECTION 16: Other information**Abbreviations and acronyms**

Skin Irrit: Skin irritation

Eye Irrit: Eye irritation

Relevant H and EUH statements (number and full text)

H315 Causes skin irritation.

H319 Causes serious eye irritation.