

Safety Data Sheet

according to UK REACH Regulation

(S)-Rolipram

Revision date: 31.05.2022

Product code: 6010

Page 1 of 5

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

(S)-Rolipram

Further trade names

2-Pyrrolidinone, 4-[3-(cyclopentyloxy)-4-methoxyphenyl]-, (4S)-

Synonyms:

(4S)-4-[3-(Cyclopentoxo)-4-methoxyphenyl]pyrrolidin-2-one

CAS No: 85416-73-5

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**

Reference standard for [11C]-(S)-Rolipram

1.3. Details of the supplier of the safety data sheet

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
e-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
e-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****GB CLP Regulation**

Acute Tox. 4; H312

Acute Tox. 4; H302

Full text of hazard statements: see SECTION 16.

2.2. Label elements**GB CLP Regulation****Hazard components for labelling**

(S)-Rolipram

Signal word: Warning**Pictograms:****Hazard statements**

H302+H312 Harmful if swallowed or in contact with skin.

Precautionary statements

P270	Do not eat, drink or smoke when using this product.
P280	Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

Safety Data Sheet

according to UK REACH Regulation

(S)-Rolipram

Revision date: 31.05.2022

Product code: 6010

Page 2 of 5

P301+P312 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
 P302+P352 IF ON SKIN: Wash with plenty of water.
 P312 Call a POISON CENTER/doctor if you feel unwell.
 P363 Wash contaminated clothing before reuse.

SECTION 3: Composition/information on ingredients**3.1. Substances****Chemical characterization**

2-Pyrrolidinone, 4-[3-(cyclopentyloxy)-4-hydroxyphenyl]-

Sum formula: C₁₆H₂₁NO₃

Molecular weight: 275.35

Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (GB CLP Regulation)			
85416-73-5	(S)-Rolipram			95 - < 100 %
	Acute Tox. 4, Acute Tox. 4; H312 H302			

Full text of H and EUH statements: see section 16.

Specific Conc. Limits, M-factors and ATE

CAS No	EC No	Chemical name	Quantity
		Specific Conc. Limits, M-factors and ATE	
85416-73-5		(S)-Rolipram	95 - < 100 %
		dermal: ATE = 1100 mg/kg; oral: ATE = 500 mg/kg	

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

Rinse mouth immediately and drink plenty of water. Medical treatment necessary.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**Water spray. Foam. Extinguishing powder. Carbon dioxide (CO₂).**5.2. Special hazards arising from the substance or mixture**In case of fire may be liberated: Carbon dioxide (CO₂). Carbon monoxide Nitrogen oxides (NO_x).**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

Additional information

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

Safety Data Sheet

according to UK REACH Regulation

(S)-Rolipram

Revision date: 31.05.2022

Product code: 6010

Page 3 of 5

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Wear personal protection equipment. Provide adequate ventilation. Avoid dust formation.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

Do not breathe dust. If handled uncovered, arrangements with local exhaust ventilation have to be used.

Advice on general occupational hygiene

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

Protect skin by using skin protective cream.

Further information on handling

When using do not eat, drink, smoke, sniff.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed.

Further information on storage conditions

storage temperature:

of °C: -25

up to °C: -15

Protect against: Light.

Store under (Gas): argon. Nitrogen.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

8.2. Exposure controls**Individual protection measures, such as personal protective equipment****Eye/face protection**

Eye protection: Tightly sealed safety glasses.

Hand protection

Tested protective gloves are to be worn: Single-use gloves. NBR (Nitrile rubber).

Skin protection

Body protection: Lab apron. Chemical resistant safety shoes.

Respiratory protection

With correct and proper use, and under normal conditions, breathing protection is not required.

Safety Data Sheet

according to UK REACH Regulation

(S)-Rolipram

Revision date: 31.05.2022

Product code: 6010

Page 4 of 5

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state: solid
 Colour: colourless

Changes in the physical state**Solubility in other solvents**

dimethylsulphoxide (DMSO). Methanol.

SECTION 10: Stability and reactivity**10.6. Hazardous decomposition products**Carbon dioxide. Carbon monoxide. Nitrogen oxides (NO_x).**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in GB CLP Regulation****Acute toxicity**

CAS No	Chemical name				
	Exposure route	Dose	Species	Source	Method
85416-73-5	(S)-Rolipram				
	oral	ATE 500 mg/kg			
	dermal	ATE 1100 mg/kg			

Carcinogenic/mutagenic/toxic effects for reproduction

No information available.

Practical experience

Causes serious eye irritation.

SECTION 12: Ecological information**12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Further information

Do not allow to enter into surface water or drains.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - used product

Safety Data Sheet

according to UK REACH Regulation

(S)-Rolipram

Revision date: 31.05.2022

Product code: 6010

Page 5 of 5

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded organic chemicals consisting of or containing hazardous substances; hazardous waste

List of Wastes Code - contaminated packaging

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by hazardous substances; hazardous waste

Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

SECTION 14: Transport information**Other applicable information**

No dangerous good in sense of these transport regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information**Relevant H and EUH statements (number and full text)**

H302	Harmful if swallowed.
H302+H312	Harmful if swallowed or in contact with skin.
H312	Harmful in contact with skin.