

Safety Data Sheet

according to UK REACH Regulation

BHTMEDA

Revision date: 31.05.2022

Product code: 6174

Page 1 of 4

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

BHTMEDA

Further trade names

1,2-Ethanediaminium, N1,N2-bis(2-hydroxyethyl)-N1,N1,N2,N2-tetramethyl-, chloride (1:2);

1,2-Ethanediaminium, N,N'-bis(2-hydroxyethyl)-N,N,N',N'-tetramethyl-, dichloride

Synonyms:

N,N'-Bis-(2-hydroxyethyl)-N,N,N',N'-tetra-methylethandiammonium dichlorid

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**Reference standard for byproduct of [¹⁸F]Fluoroethylcholine synthesis**1.3. Details of the supplier of the safety data sheet**

Company name:	ABX advanced biochemical compounds	
	Biomedizinische Forschungsreagenzien GmbH	
Street:	Heinrich-Gläser-Straße 10-14	
Place:	01454 Radeberg	
Telephone:	+49 3528 4041 60	Telefax: +49 3528 4041 65
e-mail:	info@abx.de	
Contact person:	Dr. Christoph Meyer	Telephone: +49 3528 4041 8732
e-mail:	meyer@abx.de	
Internet:	http://www.abx.de	

1.4. Emergency telephone number:

+49 3528 4041 60

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****GB CLP Regulation**

This substance is not classified as hazardous in accordance with GB CLP Regulation.

2.2. Label elements**GB CLP Regulation****Precautionary statements**

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

2.3. Other hazards

Warning - substance not yet tested completely.

SECTION 3: Composition/information on ingredients**3.1. Substances****Chemical characterization**

1,2-Ethanediaminium, N1,N2-bis(2-hydroxyethyl)-N1,N1,N2,N2-tetramethyl-, chloride (1:2);

1,2-Ethanediaminium, N,N'-bis(2-hydroxyethyl)-N,N,N',N'-tetramethyl-, dichloride

Sum formula: C₁₀H₂₆Cl₂N₂O₂

Molecular weight: 277.23

Hazardous components

none (according to UK REACH Regulation)

Safety Data Sheet

according to UK REACH Regulation

BHTMEDA

Revision date: 31.05.2022

Product code: 6174

Page 2 of 4

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air.

After contact with skin

After contact with skin, wash immediately with: Water. Change contaminated clothing.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

After ingestion

If swallowed, immediately drink: Water. Call a physician immediately.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**Co-ordinate fire-fighting measures to the fire surroundings. Suitable extinguishing media: Foam. Extinguishing powder. Carbon dioxide (CO₂). Atomized water.**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

Additional information

Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

SECTION 6: Accidental release measures**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up**Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

Do not breathe dust. Conditions to avoid: skin contact. Eye contact. Wear personal protection equipment. When using do not eat, drink, smoke, sniff.

Advice on general occupational hygiene

Change contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep container dry. Keep container tightly closed.

Further information on storage conditions

Recommended storage temperature:

of °C: -25

up to °C: -15

Protect against: Light.

Store under (Gas): argon. Nitrogen.

SECTION 8: Exposure controls/personal protection

Safety Data Sheet

according to UK REACH Regulation

BHTMEDA

Revision date: 31.05.2022

Product code: 6174

Page 3 of 4

8.1. Control parameters**Additional advice on limit values**

To date, no national critical limit values exist.

8.2. Exposure controls**Appropriate engineering controls**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state: solid
Colour: white

Changes in the physical state

Water solubility: easily soluble.

Solubility in other solvents

Methanol. Methanol.

SECTION 10: Stability and reactivity**10.4. Conditions to avoid**

Light.

10.6. Hazardous decomposition productsCarbon dioxide. Carbon monoxide. Nitrogen oxides (NO_x). Fluorines-carbon-hydrogens. Sulfur oxides.**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in GB CLP Regulation****Information on likely routes of exposure**

Toxicological data are not available.

Further information

Toxicological data are not available.

SECTION 12: Ecological information**12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

Further information

Do not allow to enter into surface water or drains.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

Contaminated packaging

Non-contaminated packages may be recycled.

SECTION 14: Transport information

Safety Data Sheet

according to UK REACH Regulation

BHTMEDA

Revision date: 31.05.2022

Product code: 6174

Page 4 of 4

Other applicable information

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Water hazard class (D): 3 - highly hazardous to water

SECTION 16: Other information