

# Safety Data Sheet

according to UK REACH Regulation

## Acetonitrile for TSTU

Revision date: 11.11.2022

Product code: PEFB-0029-V4

Page 1 of 6

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Acetonitrile for TSTU

#### Further trade names

Synonyme:

Methyl cyanide; ACN

CAS No: 75-05-8  
Index No: 608-001-00-3  
EC No: 200-835-2

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.3. Details of the supplier of the safety data sheet

Company name: ABX advanced biochemical compounds  
Biomedizinische Forschungsreagenzien GmbH  
Street: Heinrich-Gläser-Straße 10-14  
Place: 01454 Radeberg  
Telephone: +49 3528 4041 60  
e-mail: info@abx.de  
Contact person: Dr. Christoph Meyer  
e-mail: meyer@abx.de  
Internet: http://www.abx.de  
Telefax: +49 3528 4041 65  
Telephone: +49 3528 4041 8732

#### 1.4. Emergency telephone number:

+49 3528 4041 60

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### GB CLP Regulation

Flam. Liq. 2; H225  
Acute Tox. 4; H332  
Acute Tox. 4; H312  
Acute Tox. 4; H302  
Eye Irrit. 2; H319

Full text of hazard statements: see SECTION 16.

#### 2.2. Label elements

##### GB CLP Regulation

Signal word: Danger

Pictograms:



##### Hazard statements

H225 Highly flammable liquid and vapour.  
H302+H312+H332 Harmful if swallowed, in contact with skin or if inhaled.  
H319 Causes serious eye irritation.

##### Precautionary statements

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
P233 Keep container tightly closed.

**Safety Data Sheet**

according to UK REACH Regulation

**Acetonitrile for TSTU**

Revision date: 11.11.2022

Product code: PEFB-0029-V4

Page 2 of 6

P240	Ground and bond container and receiving equipment.
P241	Use explosion-proof electrical/ventilating/lighting equipment.
P242	Use non-sparking tools.
P243	Take action to prevent static discharges.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P270	Do not eat, drink or smoke when using this product.
P271	Use only outdoors or in a well-ventilated area.
P280	Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
P301+P312	IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
P330	Rinse mouth.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
P312	Call a POISON CENTER/doctor if you feel unwell.
P362+P364	Take off contaminated clothing and wash it before reuse.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P312	Call a POISON CENTER/doctor if you feel unwell.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.
P403+P235	Store in a well-ventilated place. Keep cool.

**SECTION 3: Composition/information on ingredients****3.1. Substances**Sum formula: C<sub>2</sub>H<sub>3</sub>N**Hazardous components**

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (GB CLP Regulation)			
75-05-8	acetonitrile; cyanomethane			100 %
	200-835-2	608-001-00-3		
	Flam. Liq. 2, Acute Tox. 4, Acute Tox. 4, Acute Tox. 4, Eye Irrit. 2; H225 H332 H312 H302 H319			

Full text of H and EUH statements: see section 16.

**Specific Conc. Limits, M-factors and ATE**

CAS No	EC No	Chemical name	Quantity
	Specific Conc. Limits, M-factors and ATE		
75-05-8	200-835-2	acetonitrile; cyanomethane	100 %
	inhalation: LC50 = 11 mg/l (vapours); inhalation: LC50 = 11 mg/l (dusts or mists); dermal: LD50 = 988 mg/kg; oral: ATE = 500 mg/kg		

**SECTION 4: First aid measures****4.1. Description of first aid measures****After inhalation**

Provide fresh air. In all cases of doubt, or when symptoms persist, seek medical advice.

**After contact with skin**

After contact with skin, wash immediately with: Water. Change contaminated clothing. Medical treatment necessary.

**After contact with eyes**

If product gets into the eye, keep eyelid open and rinse immediately with large quantities of water, for at least 5 minutes. Subsequently consult an ophthalmologist.

# Safety Data Sheet

according to UK REACH Regulation

## Acetonitrile for TSTU

Revision date: 11.11.2022

Product code: PEFB-0029-V4

Page 3 of 6

### After ingestion

Rinse mouth immediately and drink plenty of water. Medical treatment necessary.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable extinguishing media

Water. Carbon dioxide (CO<sub>2</sub>). Foam. Extinguishing powder.

### 5.2. Special hazards arising from the substance or mixture

Combustible. Vapours may form explosive mixtures with air.

### Additional information

Use water spray jet to protect personnel and to cool endangered containers. Suppress gases/vapours/mists with water spray jet. Contaminated fire-fighting water must be collected separately. Do not allow to enter into surface water or drains.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### General advice

Remove all sources of ignition. Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Wear personal protection equipment.

### 6.2. Environmental precautions

Do not allow to enter into surface water or drains. Explosion hazard.

### 6.3. Methods and material for containment and cleaning up

#### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advice on safe handling

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Advice on protection against fire and explosion

Keep away from sources of ignition - No smoking.

#### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Protect skin by using skin protective cream. After work, wash hands and face. When using do not eat or drink.

### 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Keep container in a well-ventilated place.

#### Hints on joint storage

Do not store together with: Material, rich in oxygen, oxidizing.

#### Further information on storage conditions

Storage at below °C: +30

Protect against: Light.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

**Safety Data Sheet**

according to UK REACH Regulation

**Acetonitrile for TSTU**

Revision date: 11.11.2022

Product code: PEFB-0029-V4

Page 4 of 6

**Exposure limits (EH40)**

CAS No	Substance	ppm	mg/m <sup>3</sup>	fibres/ml	Category	Origin
75-05-8	Acetonitrile	40	68		TWA (8 h)	WEL
		60	102		STEL (15 min)	WEL

**8.2. Exposure controls****Appropriate engineering controls**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state:	liquid
Colour:	colourless
<b>Changes in the physical state</b>	
Melting point/freezing point:	-45,7 °C
Boiling point or initial boiling point and boiling range:	81,6 °C
Flash point:	2 °C
Lower explosion limits:	3 vol. %
Upper explosion limits:	17 vol. %
Viscosity / dynamic: (at 25 °C)	0.316 mPa·s
Water solubility: (at 20 °C)	easily soluble.
Vapour pressure: (at 20 °C)	97 hPa
Density (at 20 °C):	0,786 g/cm <sup>3</sup>

**SECTION 10: Stability and reactivity****10.4. Conditions to avoid**

Keep away from heat. Ignition hazard.

**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in GB CLP Regulation****Acute toxicity**

Acute toxicity, oral. Acute toxicity, inhalant. Acute toxicity, dermal.

**Safety Data Sheet**

according to UK REACH Regulation

**Acetonitrile for TSTU**

Revision date: 11.11.2022

Product code: PEFB-0029-V4

Page 5 of 6

CAS No	Chemical name				
	Exposure route	Dose	Species	Source	Method
75-05-8	acetonitrile; cyanomethane				
	oral	ATE mg/kg	500		
	dermal	LD50 mg/kg	988	Rabbit	IUCLID
	inhalation vapour	LC50	11 mg/l		
	inhalation dust/mist	LC50	11 mg/l		

**Irritation and corrosivity**

Irritating to eyes.

**SECTION 12: Ecological information****12.1. Toxicity**

CAS No	Chemical name					
	Aquatic toxicity	Dose	[h]   [d]	Species	Source	Method
75-05-8	acetonitrile; cyanomethane					
	Acute fish toxicity	LC50 mg/l	1640	96 h	Pimephales promelas	IUCLID

**12.3. Bioaccumulative potential****Partition coefficient n-octanol/water**

CAS No	Chemical name	Log Pow
75-05-8	acetonitrile; cyanomethane	-0,34

**12.5. Results of PBT and vPvB assessment**

This substance does not meet the PBT/vPvB criteria of UK REACH.

**12.6. Endocrine disrupting properties**

This substance does not have endocrine disrupting properties with respect to non-target organisms.

**Further information**

Do not allow to enter into surface water or drains.

**SECTION 13: Disposal considerations****13.1. Waste treatment methods****Disposal recommendations**

Dispose of waste according to applicable legislation.

**Contaminated packaging**

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

**SECTION 14: Transport information****Land transport (ADR/RID)**

<b>14.1. UN number or ID number:</b>	UN 1648
<b>14.2. UN proper shipping name:</b>	ACETONITRILE
<b>14.3. Transport hazard class(es):</b>	3
<b>14.4. Packing group:</b>	II
Hazard label:	3

**Safety Data Sheet**

according to UK REACH Regulation

**Acetonitrile for TSTU**

Revision date: 11.11.2022

Product code: PEFB-0029-V4

Page 6 of 6

**Inland waterways transport (ADN)**

<b><u>14.1. UN number or ID number:</u></b>	UN 1648
<b><u>14.2. UN proper shipping name:</u></b>	ACETONITRILE
<b><u>14.3. Transport hazard class(es):</u></b>	3
<b><u>14.4. Packing group:</u></b>	II
Hazard label:	3

**Marine transport (IMDG)**

<b><u>14.1. UN number or ID number:</u></b>	UN 1648
<b><u>14.2. UN proper shipping name:</u></b>	ACETONITRILE
<b><u>14.3. Transport hazard class(es):</u></b>	3
<b><u>14.4. Packing group:</u></b>	II
Hazard label:	3

**Air transport (ICAO-TI/IATA-DGR)**

<b><u>14.1. UN number or ID number:</u></b>	UN 1648
<b><u>14.2. UN proper shipping name:</u></b>	ACETONITRILE
<b><u>14.3. Transport hazard class(es):</u></b>	3
<b><u>14.4. Packing group:</u></b>	II
Hazard label:	3

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information**

Restrictions on use (REACH, annex XVII):

Entry 3, Entry 40, Entry 75

2004/42/EC (VOC): 100 % (786 g/l)

**National regulatory information**

Water hazard class (D): 2 - obviously hazardous to water

**SECTION 16: Other information****Relevant H and EUH statements (number and full text)**

H225	Highly flammable liquid and vapour.
H302	Harmful if swallowed.
H302+H312+H332	Harmful if swallowed, in contact with skin or if inhaled.
H312	Harmful in contact with skin.
H319	Causes serious eye irritation.
H332	Harmful if inhaled.